

GLOBACLEAN[®]

CLEANING AND DISINFECTION AGENT FOR INSTRUMENTS

concentrate

Globaclean[®] is a cost-effective agent, excellent for removing organic dirt, with no aldehydes, phenols or chlorine, for washing and disinfection of instruments prior to high-level disinfection and sterilisation. Suitable for use in ultrasonic baths too. When used in washing machines, the adding of foam blocker is recommended.

Also conforms to requirements in food industry, beauty services etc.

PRODUCT PROPERTIES

- wide spectrum of efficacy
- suitable for use in instrument and ultrasonic baths as well as endoscope washing machines (recommended to add foam blocker)
- not tuberculocidal
- ready-to-use working solution can be preserved in a closed vessel for 4 weeks
- the working solution in use can be used for 7 days

USAGE

Make the working solution of clean water and concentrate according to the exposure time needed (look at the table). Disassemble the instruments and submerge them into a lid-equipped tub filled with the working solution. Make sure all the instruments are covered with the solution. Once the exposure time has passed, wash or rinse the instruments under running water (at least tap water quality) or with sterile water, as needed and dry with lint-free towel or disposable tissues (instruments can't be left out to dry in the air). It's not recommended to keep the instruments/accessories in the solution longer than the exposure time.

In order to make a working solution out of Globaclean[®], use water that is cool or at room temperature.

Note! Do not immerse dental handpieces and turbines in the solution. Process them with a rapid disinfectant (like Vitasept[®] Wipes med) by wiping their outer surfaces.

WORKING SOLUTION CONCENTRATIONS AND EXPOSURE

Less soiled instruments	1%	10ml/l	5min
Grossly soiled instruments	3%	30ml/l	5min
Grossly soiled instruments (e.g HBV)	3%	30ml/l	30min

MICROBIOLOGICAL SPECTRUM

Product has bactericidal, fungicidal and virucidal (incl. HBV, HIV, HCV, coronavirus, influenza, vaccinia virus) efficacy.

CONFORMITY TO STANDARDS

EN 1040, EN 1276, EN 13624, EN 13727, EN 14476, EN 14561, EN 14562, DVG, DVV/RKI, DGHM/ÖGHMP, AFNOR NFT

STORAGE

See the best before date on the package. Store at temperatures 5-25°C, keep away from direct sunlight, avoid freezing to assure the quality of the product.

UTILIZATION

The working solution is ecologically safe, all active ingredients are biodegradable in flush water. Plastic packaging materials are recyclable. For collection, the plastic packaging is suitable when rinsed, the bottle and cap separated. The packaging box is sorted for cardboard collection. The compatibility of the packaging and the product is quality assured.

SAFETY INFORMATION

See the product safety data sheet for additional safety information.

COMPOSITION

100g of the product contains active ingredients: didecyldimethylammonium chloride 4,5g

pH 12


Product bears the CE marking according to Medical Device directive 93/42/EEC



PACKAGES



1L dosing bottle



5L canister



Goodpoint Chemicals OÜ
Urda tee 2/1, Jälgimäe 76404,
Saku vald, Harjumaa, Estonia
www.goodpointchemicals.com
info@goodpointchemicals.com

